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- Elias JM, Gown AM, Nakamura RM, Wilbur DC, Herman GE, Jaffe ES, Battifora H, Brigati DJ. Quality control in immunohistochemistry. Report of a workshop sponsored by the Biological Stain Commission. Am J Clin Pathol 92 (6): 836-843, 1989.
- Sheehan DC, Hrapchak BB. Theory and practice of histotechnology, 2nd Edition. The C.V. Mosby Company, St. Louis, 1980.
- College of American Pathologists Laboratory Accreditation Program, Anatomic Pathology Checklist, 2001.
- CLSI. Quality Assurance for Immunocytochemistry: Approved Guideline. CLSI document MM4-A- (ISBN 1-56238-396-5). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 1999.
- Herman GE, Elfont EA. The taming of immunohistochemistry: the new era of quality control. Biotech Histochem 66(4): 194-199, 1991.
- Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immunohistochemistry. Am J Clin Pathol 73(5): 626-32, 1980.
- Nadji M, Morales AR. Immunoperoxidase: part 1. The technique and its pitfalls. Lab Med 14: 767, 1983.
- Roche PC, Hsi ED. Immunohistochemistry-Principles and Advances. Manual of Clinical Laboratory Immunology, 6th edition. (NR Rose Ed.) ASM Press, 2002.
- Taylor C, Cote RJ. Immunomicroscopy: A Diagnostic Tool for the Surgical Pathologist, 2nd Edition. W.B. Saunders Company, Philadelphia. 1986.

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CE

16

ultraView™ Universal DAB Detection Kit

Catalog number 760-500

INDICATIONS AND USE

Intended Use

This detection kit is intended for *in vitro* diagnostic use.

Ventana Medical Systems' (Ventana) *ultraView* Universal DAB Detection Kit is an indirect, biotin-free system for detecting mouse IgG, mouse IgM and rabbit primary antibodies. The kit is intended to identify targets by immunohistochemistry (IHC) in sections of formalin fixed, paraffin embedded and frozen tissue that are stained on the Ventana NexES® and BenchMark® Series automated slide stainers. The clinical interpretation of any staining, or the absence of staining, must be complemented by morphological studies and evaluation of proper controls. Evaluation must be made by a qualified pathologist within the context of the patient's clinical history and other diagnostic tests.

Summary and Explanation

Immunohistochemistry (IHC) is a technique used in laboratories for diagnostic purposes. The fundamental concept of IHC is the localization of antigens in tissue sections using specific primary antibodies. Once the antibody-antigen binding occurs it must be visualized with a colored reaction visible by light microscopy or through the use of fluorescence microscopy. The specific primary antibody may carry the signal generating molecule providing direct visualization of the binding. Alternatively, indirect methods of visualization require additional steps to localize the specific antibody and generate a signal. The most common techniques of indirect methods use a secondary antibody directed against the species of primary antibody and an enzyme with a corresponding substrate chromogen system. This combination results in a colored precipitate at the site of specific antibody binding. The *ultraView* Universal DAB Detection kit uses an indirect method to visualize specific antibodies bound to antigens by depositing a brown colored precipitate.

Principles and Procedures

The *ultraView* Universal DAB Detection Kit detects specific mouse and rabbit primary antibodies bound to an antigen in paraffin embedded or frozen tissue sections. The specific antibody is located by an enzyme labeled secondary antibody. The complex is then visualized utilizing a precipitating enzyme product.

The *ultraView* Universal DAB Detection Kit utilizes a cocktail of enzyme labeled secondary antibodies that locate the bound primary antibody. The complex is then visualized with hydrogen peroxide substrate and 3, 3' - diaminobenzidine tetrahydrochloride (DAB) chromogen, which produce a dark brown precipitate that is readily detected by light microscopy. Each step is incubated for a precise time and temperature. At the end of each incubation step, the Ventana automated slide stainer washes the sections to remove unbound material and applies a liquid coverslip which minimizes the evaporation of the aqueous reagents from the slide.¹ Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

For more detailed information on instrument operation, refer to the appropriate Ventana automated slide stainer Operator's Manual.

MATERIALS AND METHODS

Reagents Provided

ultraView Universal DAB Detection Kit contains sufficient reagent for 250 tests and contains:

- 1 – 25 ml dispenser Inhibitor; contains 3% hydrogen peroxide solution
- 1 – 25 ml dispenser Universal HRP Multimer; contains a cocktail of HRP labeled antibodies (goat anti-mouse IgG, goat anti-mouse IgM, and goat anti-rabbit) (<50 µg/mL) in a buffer containing protein with ProClin® 300 preservative*
- 1 – 25 ml dispenser DAB Chromogen; contains 3, 3' - diaminobenzidine tetrahydrochloride in a proprietary stabilizer solution with a proprietary preservative
- 1 – 25 ml dispenser DAB H₂O₂; contains 0.04% hydrogen peroxide in a phosphate buffer solution
- 1 – 25 ml dispenser Copper; contains copper sulfate (5g/l) in an acetate buffer with proprietary preservative

*For preservative information please refer to the "Warnings and Precautions" section.

Reconstitution, Mixing, Dilution, Titration

The detection kit is optimized for use on a Ventana automated slide stainer. No reconstitution, mixing, dilution, or titration of kit reagents is required. Further dilution may result in loss of antigen staining. The user must verify any such changes. Differences in tissue processing and technical procedure in the user's laboratory may produce significant variability in results and require regular use of controls (see Quality Control Procedures section).

Materials and Reagents Needed But Not Provided

The following reagents and materials may be required for staining but are not provided:

1. Negative reagent control (Ventana Negative Control Mouse Ig (Cat No 760-2014), CONFIRM™ Negative Control Rabbit Ig (Cat No 760-1029)
2. Superfrost™ Plus (VWR Cat. No. 48311-703 or equivalent) microscope slides
3. Positive and negative tissue controls (consult antibody package insert for recommended types)
4. Drying oven capable of maintaining a temperature of 60° C ± 5° C
5. Bar code labels (appropriate for negative reagent control and primary antibody being tested)
6. Staining jars or baths
7. Timer (capable of 2-10 minute intervals)
8. Xylene (Histological grade)
9. Ethanol or reagent alcohol (Histological grade)
 - 100% solution: undiluted ethanol or reagent alcohol
 - 95% solution: mix 95 parts ethanol or reagent alcohol with 5 parts of deionized water
 - 80% solution: mix 80 parts ethanol or reagent alcohol with 20 parts deionized water
10. Deionized or distilled water
11. Biocare Medical's Decloaking Chamber (Cat No DC2002) (NexES IHC automated slide stainers)*
12. Tissue-Tek® staining dishes*
13. NexES IHC, BenchMark Series automated slide stainers
14. Primary antibody
15. Ventana Antibody Diluent* (Cat No 251-018)
16. Ventana APK Wash (10X) (Cat. No. 250-042) (NexES IHC automated slide stainers)
17. Ventana Liquid Coverslip™ (Low Temperature) (Cat. No. 250-009) (NexES IHC automated slide stainers)
18. Ventana EZ Prep™ (10X) (Cat. No. 950-102) (BenchMark Series automated slide stainers)
19. Ventana Reaction Buffer (10X) (Cat. No. 950-300) (BenchMark Series automated slide stainers)
20. Ventana Liquid Coverslip (High Temperature) (Cat. No. 650-010) (BenchMark Series automated slide stainers)
21. Ventana Cell Conditioning 1 (Pre-dilute) (Cat. No. 950-124) or Cell Conditioning 2 (Pre-dilute) (Cat No 950-123) (BenchMark Series automated slide stainers)*
22. Ventana Protease 1 (Cat No 760-2018), Protease 2 (Cat No 760-2019), or Protease 3 (Cat No 760-2020)*
23. Ventana Hematoxylin counterstain* (Cat No 760-2021)
24. Ventana Hematoxylin II counterstain* (Cat. No. 790-2208)
25. Ventana Bluing Reagent* (Cat. No. 760-2037)
26. Permanent Mounting Medium (Permount®, Fisher Cat. No. SP15-500 or equivalent)
27. Cover glass (sufficient to cover tissue such as VWR Cat. No. 48393-60 or equivalent)
28. Light microscope (20-80X)
29. Absorbent wipes*
30. Automated coverslipper (such as Tissue-Tek® SCA automated coverslipper)
 - *As needed for specific applications.

Storage and Handling

Store at 2° C to 8° C. Do not freeze. The user must validate any storage conditions other than those specified in the package insert. This detection kit can be used immediately after removal from the refrigerator.

To ensure proper reagent delivery and stability of each reagent, after every run the cap must be replaced and the dispenser must be immediately placed in the refrigerator in an upright position.

Every detection kit is expiration dated. When properly stored, the reagents are stable to the date indicated on the label. Do not use detection kit beyond the expiration date.

There are no definitive signs to indicate instability of this product; therefore, positive and negative controls should be run simultaneously with unknown specimens. Your local Ventana office should be contacted immediately if there is an indication of reagent instability.

Specimen Collection and Preparation for Analysis

Formalin fixed, paraffin embedded tissues and frozen tissues are suitable for use with *ultraView* Universal DAB Detection Kit and a Ventana automated slide stainer (see Materials and Reagents Needed, But Not Provided section). The recommended tissue fixative is 10% neutral buffered formalin.² Variable results may occur as a result of prolonged fixation or special processes such as decalcification of bone marrow preparations.

Each section should be cut to the appropriate thickness (2-5µ) for the primary antibody being used and placed on a glass slide. Slides containing the tissue section may be baked/heated for at least 2 hours (but not longer than 24 hours) in a 60° C ± 5° C oven. Slide heating is used to dry the tissue post slide mounting and to enhance tissue adhesion to the glass. Extended heating of the tissue might result in decreased antigen availability. Consult the primary antibody package insert to identify heating limitations.

Properly fixed and embedded tissues expressing the antigen will remain stable for at least 2 years if stored in a cool location (15° C - 25° C). The Clinical Laboratory Improvement Act (CLIA) of 1988, 42CFR493.1259 (b) requires that "The laboratory must retain slides at least

ten years from the date of examination and retain specimen blocks at least two years from date of examination".
Frozen sections should be cut to the appropriate thickness (4-5µ) for the primary antibody being used, picked up on a glass microscope slide and immediately placed in cold acetone (4° C to 8° C) for ten minutes. The sections are then air dried for a minimum of thirty minutes and preferably overnight. Apply appropriate barcode label to dry slides.

Manual Deparaffinization Procedure

Required when using the NexES IHC automated slide stainer or if deparaffinization is not selected on the BenchMark Series automated slide stainer:

1. For instructions on when to label slides with barcode label, refer to the Instructions for Use section of the specific automated slide stainer Operator's Manual.
2. Immerse the slides sequentially in 3 xylene baths for 5 ± 1 minutes each.
3. Transfer the slides to 100% ethanol and immerse sequentially in 2 baths for 3 ± 1 minutes each.
4. Transfer the slides to 95% ethanol and immerse them in a bath of this solution for 3 ± 1 minutes.
5. Transfer the slides to 80% ethanol and immerse them in this solution for 3 ± 1 minutes.
6. Transfer the slides to a bath of deionized or distilled water and dip a minimum of 10 times.
7. Transfer slides to APK Wash (1X) or buffer solution as appropriate. For APK Wash, the slides should remain in this solution until you are ready to perform the staining run. For buffer solution, the slides should remain in this solution until you are ready to perform the antigen unmasking procedure. Do not allow the slides to dry.

Slides stained on the BenchMark Series automated slide stainers can be deparaffinized on the instrument. If this option is selected, apply barcode labels to slides and place slides on the instrument. If the option is not selected, follow the Manual Deparaffinization Procedure above.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use:

1. Take reasonable precautions when handling reagents. Use disposable gloves when handling suspected carcinogens or toxic materials (example: xylene or formaldehyde). Do not use near open flame.
2. Avoid contact of reagents with eyes and mucous membranes. If reagents come in contact with sensitive areas, wash with copious amounts of water.
3. Patient specimens and all materials contacting them should be handled as biohazardous materials and disposed of with proper precautions. Never pipette by mouth.
4. Avoid microbial contamination of reagents, as this could produce incorrect results.
5. Incubation times and temperatures other than those specified may give erroneous results. The user must validate any such change.
6. The reagents have been optimally diluted, and further dilution may result in loss of antigen staining. The user must validate any such change. Do not mix reagent dispensers from multiple lots of the product.
7. Symptoms of overexposure to ProClin® 300 include skin and eye irritation, and irritation of mucous membranes and upper respiratory tract. The concentration of ProClin® 300 in this product is less than or equal to 0.05% and does not meet the OSHA criteria for a hazardous substance. Systemic allergic reactions are possible in sensitive individuals. Proprietary preservatives and stabilizers do not meet the OSHA criteria for a hazardous substance.
8. Possible carcinogen. The National Toxicology Program has listed benzidine, a closely related compound to 3, 3' - diaminobenzidine tetrahydrochloride (DAB), as a known human carcinogen.
9. Concentrated forms of propylene glycol have been associated with teratogenic effects in laboratory animals. Use disposable gloves and take reasonable precautions when handling.
10. Consult local or state authorities with regard to recommended method of disposal.

INSTRUCTIONS FOR USE

Step by Step Procedure

Ventana *ultraView* Universal DAB Detection Kit has been developed for use on Ventana automated slide stainers in combination with Ventana primary antibodies and accessories. The parameters for the automated procedures can be displayed, printed and edited according to the procedure in the Operator's Manual for the individual automated slide stainer. Other operating parameters for the automated slide stainers have been preset at the factory.

The procedures for staining on the Ventana automated slide stainers are as follows. For more detailed instructions and additional protocol options refer to your Operator's Manual. Whether a sample requires Antigen Unmasking is antibody dependent. Please check the antibody package insert for directions.

NexES IHC Automated Slide Stainers

If Manual Antigen Unmasking Is Required:

1. Slides are to be deparaffinized through a series of xylene and gradient alcohols to water and then to appropriate buffer (as indicated by the antibody package insert).

Perform antigen unmasking procedure (as indicated by the antibody package insert) and transfer slides to APK Wash (1X).

2. Load the primary antibody, appropriate detection kit dispensers and required accessory reagents onto the reagent tray and place the reagent tray on the automated slide stainer. Check bulk fluids and waste.
3. Dry the painted end of the slide and then apply slide barcode label that corresponds to the antibody protocol to be performed.
4. Load the deparaffinized, antigen unmasked, labeled slides from the APK Wash (1X) onto the NexES IHC Automated Slide Stainer. Slides should remain in this solution until ready to perform staining run.

NexES IHC Automated Slide Stainers and Manual

If Antigen Unmasking Is Not Required:

1. Apply barcode labels to the slide. Slides are to be deparaffinized through a series of xylene and gradient alcohols to water and then to APK Wash (1X).
2. Load the primary antibody, appropriate detection kit dispensers and required accessory reagents onto the reagent tray and place the reagent tray on the automated slide stainer. Check bulk fluids and waste.
3. Load the deparaffinized, labeled slides from the APK Wash (1X) onto the NexES IHC Automated Slide Stainer. Slides should remain in this solution until ready to perform staining run.

BenchMark Series Automated Slide Stainers

1. Apply slide barcode label that corresponds to the primary antibody protocol to be performed.
2. Load the primary antibody, appropriate detection kit dispensers and required accessory reagents onto the reagent tray and place the reagent tray on the automated slide stainer. Check bulk fluids and waste.
3. Load the slides onto the automated slide stainer.

For All Instruments

1. Start the staining run.
2. At the completion of the run, remove the slides from the automated slide stainer.
3. Wash slides in an aqueous solution containing a mild dishwashing detergent or dilute alcohol to remove the coverslip solution; dehydrate and coverslip with permanent mounting media in the usual manner.
4. Slides may be read immediately upon removal from the staining platform. Stained slides are stable for at least two years, if properly stored at room temperature (15° C to 25° C).

Quality Control Procedures

Positive Tissue Control

A positive tissue control must be run with every staining procedure performed. The positive staining tissue components are used to confirm that the antibody was applied and the instrument functioned properly. This tissue may contain both positive and negative staining cells or tissue components and may serve as both the positive and negative control tissue. Control tissues should be fresh autopsy, biopsy or surgical specimens prepared or fixed as soon as possible in a manner identical to the test sections. Such tissues may monitor all steps of the procedure, from tissue preparation through staining. Use of a tissue section fixed or processed differently from the test specimen will provide control for all reagents and method steps except fixation and tissue processing.

A tissue with weak positive staining is more suitable for optimal quality control and for detecting minor levels of reagent degradation.

Known positive tissue controls should be utilized only for monitoring the correct performance of processed tissues and test reagents, not as an aid in determining a specific diagnosis of patient samples. If the positive tissue controls fail to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control

The same tissue used for the positive tissue control may be used as the negative tissue control. The variety of cell types present in most tissue sections offers internal negative control sites, but this must be verified by the user. The components that do not stain should demonstrate the absence of specific staining, and provide an indication of background staining. If specific staining occurs in the negative tissue control sites, results with the patient specimens should be considered invalid.

Unexplained Discrepancies

Unexplained discrepancies in controls should be referred to your local Ventana office immediately. If quality control results do not meet specifications, patient results are invalid. If discrepancies occur, refer to the Troubleshooting section of this insert. Identify and correct the problem, then repeat the patient samples.

Negative Reagent Control

A negative reagent control must be run for every specimen to aid in the interpretation of results. A negative reagent control is used in place of the primary antibody to evaluate nonspecific staining. The slide should be stained with Negative Control Mouse Ig, CONFIRM Negative Control Rabbit Ig, or Rabbit Negative Control as appropriate. If an alternative negative reagent control is used, dilute to the same concentration as the primary

18

antibody antiserum with Ventana Antibody Diluent. The diluent alone may be used as an alternative to the previously described negative reagent controls. The incubation period for the negative reagent control should equal the primary antibody incubation period. When panels of several antibodies are used on serial sections, a negative reagent control on one slide may serve as a negative or nonspecific binding background control for other antibodies.

Assay Verification

Prior to initial use of a primary antibody or staining system in a diagnostic procedure, the specificity of the primary antibody should be verified by testing on a series of tissues with known immunohistochemistry performance characteristics representing known positive and negative tissues (refer to the Quality Control Procedures listed in the primary antibody product insert and to the Quality Control recommendations of the College of American Pathologists Laboratory Accreditation Program, Anatomic Pathology Checklist³, or the CLSI Approved Guideline⁴ or both documents). These quality control procedures should be repeated for each new antibody lot, or whenever there is a change in assay parameters.

Interpretation of Results

The Ventana *ultraView* Universal DAB Detection Kit causes a dark brown colored reaction product to precipitate at the antigen sites localized by the primary antibody. A qualified pathologist experienced in immunohistochemical procedures must evaluate controls and qualify the stained product before interpreting results. Staining of negative controls must be noted first, and these results compared to the stained material to verify that the signal generated is not the cause of nonspecific interactions.

Positive Tissue Control

The stained positive tissue control should be examined first to ascertain that all reagents are functioning properly. The presence of an appropriately colored reaction product within the target cells is indicative of positive reactivity. Depending on the incubation length and potency of the hematoxylin used, counterstaining will result in a pale to dark blue coloration of cell nuclei. Excessive or incomplete counterstaining may compromise proper interpretation of results.

If the positive tissue control fails to demonstrate positive staining, any results with the test specimens should be considered invalid.

Negative Tissue Control

The negative tissue control should be examined after the positive tissue control to verify the specific labeling of the target antigen by the primary antibody. The absence of specific staining in the negative tissue control confirms the lack of antibody cross reactivity to cells or cellular components. If specific staining occurs in the negative tissue control, results with the patient specimen should be considered invalid.

Nonspecific staining, if present, will have a diffuse appearance. Sporadic light staining of connective tissue may also be observed in sections from excessively formalin fixed tissues. Intact cells should be used for interpretation of staining results. Necrotic or degenerated cells often stain nonspecifically.

Patient Tissue

Patient specimens should be examined last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen in question was not detected, not that the antigen is absent in the cells or tissue assayed. If necessary, use a panel of antibodies to aid in the identification of false negative reactions. The morphology of each tissue sample should also be examined utilizing a hematoxylin and eosin stained section when interpreting any immunohistochemical result. The patient's morphologic findings and pertinent clinical data must be interpreted by a qualified pathologist.

LIMITATIONS

General Limitations

1. Immunohistochemistry is a multiple step diagnostic process that requires specialized training in the selection of the appropriate reagents and tissues, fixation, processing, preparation of the immunohistochemistry slide, and interpretation of the staining results.
2. Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be a consequence of variations in fixation and embedding methods, or inherent irregularities within the tissue.
3. Excessive or incomplete counterstaining may compromise proper interpretation of results.
4. The clinical interpretation of any positive staining, or its absence, must be evaluated within the context of clinical history, morphology and other histopathological criteria. The clinical interpretation of any staining, or its absence, must be complemented by morphological studies and proper controls as well as other diagnostic tests. It is the responsibility of a qualified pathologist to be familiar with the antibodies, reagents and methods used to interpret the stained preparation. Staining must be performed in a certified licensed laboratory under the supervision of a pathologist who is responsible

for reviewing the stained slides and assuring the adequacy of positive and negative controls.

5. Ventana provides reagents at optimal dilution for use when the provided instructions are followed. Any deviation from recommended test procedures may invalidate expected results. Appropriate controls must be employed and documented. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results.
6. Reagents may demonstrate unexpected reactions in previously untested tissues. The possibility of unexpected reactions even in tested tissue groups cannot be completely eliminated because of biological variability of antigen expression in neoplasms, or other pathological tissues.⁵ Contact your local Ventana office with documented unexpected reactions.
7. Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HBsAg) may exhibit nonspecific staining with horseradish peroxidase.⁶
8. When used in blocking steps, normal sera from the same animal source as the secondary antisera may cause false negative or false positive results due to autoantibodies or natural antibodies.
9. False positive results may be seen because of nonimmunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes), endogenous peroxidase activity (cytochrome C).⁷
10. As with any immunohistochemistry test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells or tissue assayed.

Specific Limitations

1. Each step of the detection kit procedure has been optimized on the Ventana automated slide stainers and is preset. Because of variation in tissue fixation and processing, it may be necessary to increase or decrease the primary antibody incubation time on individual specimens. For further information on fixation variables, refer to "Immunohistochemistry Principles and Advances"⁸ or "Immunomicroscopy: A Diagnostic Tool for the Surgical Pathologist".⁹
2. The detection kit, in combination with Ventana primary antibodies and accessories, detects antigen that survives routine formalin, tissue processing and sectioning. Users who deviate from recommended test procedures are responsible for interpretation and validation of patient results.

SUMMARY OF EXPECTED RESULTS

Reproducibility Studies

ultraView Universal DAB Detection Kit reproducibility testing was performed by staining serial sections from 3 neutral buffered formalin fixed, paraffin embedded tissues using 3 primary antibodies, a mouse IgG (anti-Ki67) stained on breast carcinoma, a mouse IgM (anti-CD15) stained on a Hodgkin's Lymphoma xenograft and a rabbit IgG (anti-S100) on a melanoma using BenchMark, BenchMark XT and NexES IHC automated slide stainers. All primary antibodies were incubated for 16 minutes and slides counterstained using Hematoxylin II followed by Bluing Reagent. All slides stained with a primary antibody were compared against each other for staining appropriateness and intensity and scored by a qualified slide reader.

1. Intra run reproducibility (same primary antibody stained on a platform compared) staining runs were performed one per day on 3 separate days using 3 different stainers: BenchMark, BenchMark XT, & NexES IHC for a total of 9 staining runs. Intra run reproducibility for BenchMark XT was 100% (30 of 30 slides per primary antibody per run, a total of 90 slides stained per antibody across the 3 runs) and 100% (20 of 20 slides per primary antibody per run, a total of 60 slides per antibody across 3 runs for each platform: BenchMark and NexES IHC).
2. Inter run reproducibility was calculated based on the number of slides stained in 3 runs per stainer type. Staining runs were performed one per day on 3 separate days using 3 different stainers: BenchMark, BenchMark XT, & NexES IHC. Inter run reproducibility for BenchMark XT was 100% (90 of 90 slides, 30 slides for each primary antibody across 3 separate runs) and 100% (60 of 60 slides, 20 slides for each primary antibody across 3 separate runs for each platform) for BenchMark and NexES IHC.
3. Inter instrument reproducibility was calculated based on the number of slides stained in 9 runs across all stainer types. Staining runs were performed 1 per day on 3 separate days using 3 different stainers: BenchMark, BenchMark XT, & NexES IHC. *ultraView* Universal DAB Detection Kit inter instrument reproducibility is 100% (210 of 210 stained slides, evaluated slides included all 3 primary antibodies).

The performance of the *ultraView* Universal DAB Detection Kit was also determined by staining 53 normal and disease state neutral buffer formalin fixed tissues using 20 primary antibodies. The selected antibodies represented species and isotypes commonly used to perform IHC. In 100% of the cases tested (53 of 53 cases stained) *ultraView* Universal DAB Detection Kit recognized the primary antibody and provided the appropriate clinical utility as determined by the qualified slide reader.

TROUBLESHOOTING

1. If the positive control exhibits weaker staining than expected, other positive controls run concurrently should be checked to determine if it is due to the primary antibody or one of the common secondary reagents.

2. If the positive control is negative, it should be checked to ensure that the slide has the proper barcode label. If the slide is labeled properly, other positive controls run concurrently should be checked to determine if it is due to the primary antibody or one of the common secondary reagents. Tissues may have been improperly collected, fixed or deparaffinized. The proper procedure should be followed for collection, storage and fixation.
3. If all of the paraffin has not been removed, there may be no staining. The deparaffinization procedure should be repeated.
4. If specific antibody staining is too intense, the run should be repeated with the primary antibody incubation time shortened by 4 minute intervals until the desired stain intensity is achieved.
5. If tissue sections wash off the slide, slides should be checked to ensure that they are positively charged.
6. For corrective action, refer to the Step By Step Procedure section, the automated slide stainer Operator's Manual or contact your local Ventana office.

REFERENCES

1. Elias JM, Gown AM, Nakamura RM, Wilbur DC, Herman GE, Jaffe ES, Battifora H, Brigati DJ. Quality control in immunohistochemistry. Report of a workshop sponsored by the Biological Stain Commission. Am J Clin Pathol 92 (6): 836-843, 1989.
2. Sheehan DC, Hrapchak BB. Theory and practice of histotechnology, 2nd Edition. The C.V. Mosby Company, St. Louis, 1980.
3. College of American Pathologists Laboratory Accreditation Program, Anatomic Pathology Checklist, 2001.
4. CLSI. Quality Assurance for Immunocytochemistry: Approved Guideline. CLSI document MM4-A- (ISBN 1-56238-396-5). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 1999.
5. Herman GE, Elfont EA. The taming of immunohistochemistry: the new era of quality control. Biotech Histochem 66(4): 194-199, 1991.
6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immunohistochemistry. Am J Clin Pathol 73(5): 626-32, 1980.
7. Nadji M, Morales AR. Immunoperoxidase: part 1. The technique and its pitfalls. Lab Med 14: 767, 1983.
8. Roche PC, Hsi ED. Immunohistochemistry-Principles and Advances. Manual of Clinical Laboratory Immunology, 6th edition. (NR Rose Ed.) ASM Press, 2002.
9. Taylor C, Cote RJ. Immunomicroscopy: A Diagnostic Tool for the Surgical Pathologist, 2nd Edition. W.B. Saunders Company, Philadelphia. 1986.

CONTACT INFORMATION

North America
Ventana Medical Systems, Inc.
1910 E. Innovation Park Drive
Tucson, Arizona 85755
U.S.A.
800 227 2155 (U.S. only)
+1 520 887 2155

Europe
Ventana Medical Systems, S.A.
Parc d'Innovation – BP 30144
Rue G. de Kaysersberg
F - 67404 Illkirch CEDEX
France
+33 (0) 3 90 40 52 00

EU Representative
MDCI Ltd.
Arundel House
1 Liverpool Gardens
Worthing
West Sussex BN11 1SL

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Australia, New Zealand
Ventana Medical Systems Pty Ltd
5/39 Grand Boulevard
Montmorency VIC 3094
Australia
+61 (0) 3 9431 6064

INTELLECTUAL PROPERTY

CONFIRM™, EZ Prep™, *ultraView*™, and Liquid Coverslip™ are trademarks of Ventana Medical Systems, Inc.; BenchMark®, NexES® and Ventana® are registered trademarks of Ventana Medical Systems, Inc.

ProClin® is a registered trademark of Rohm & Haas Company.

Permount® is a registered trademark of Fisher Scientific Company.

Tissue-Tek® is a registered trademark of Sakura Finetek.

Superfrost™ Plus is a trademark of Erie Scientific Company.

Ventana grants to Purchaser a single use only license under the following patents: U.S. Pat. Nos. 6,045,759; 6,192,945; 6,416,713; 6,945,128 and foreign counterparts.

ultraView Universal DAB Detection Kit covered by patents pending.

CE

20