

“巴德”威力克力育艾斯植入式注射座 “Bard” PowerPort ClearVUE isp Implantable Port

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- 空氣栓塞症 (Air Embolism)
- 過敏反應 (Allergic Reaction)
- 出血 (Bleeding)
- 臂神經叢受傷 (Brachial Plexus Injury)
- 心律不整 (Cardiac Arrhythmia)
- 心臟穿刺 (Cardiac Puncture)
- 心包填塞 (Cardiac Tamponade)
- 導管或注射座糜爛遍及皮膚 (Catheter or Port Erosion through the Skin)
- 導管栓塞 (Catheter Embolism)
- 導管阻塞 (Catheter Occlusion)
- 導管或注射座相關敗血症 (Catheter or Port related Sepsis)
- 因鎖骨與第一肋骨間的擠壓使阻塞、受損或斷裂 (Damage or Breakage due to Compression between the Clavicle and First Rib)
- 器材旋轉或擠出 (Device Rotation or Extrusion)
- 心內膜炎 (Endocarditis)
- 藥物外滲 (Extravasation)
- 纖維鞘形成 (Fibrin Sheath Formation)
- 血腫 (Hematoma)
- 血胸 (Hemothorax)
- 水胸 (Hydrothorax)
- 感染·包括但不限於注射座口·導管隧道器和/或血流 (Infection, including but not limited to pocket, catheter tunnel, and/or blood stream)
- 植入區域的皮膚發炎·壞死或留疤 (Inflammation, Necrosis, or Scarring of Skin Over Implant Area)
- 植入性器材不耐反應 (Intolerance Reaction to Implanted Device)
- 血管或內臟裂傷 (Laceration of Vessels or Viscus)
- 注射座口或周圍疼痛 (Pain at or around port pocket site)
- 血管或內臟穿孔 (Perforation of Vessels or Viscus)
- 氣胸 (Pneumothorax)
- 與局部及全身性麻醉·外科手術及術後恢復有關的風險 (Risks Normally Associated with Local or General Anesthesia, Surgery, and Post- Operative Recovery)
- 自發性導管尖端位置不正或收縮 (Spontaneous Catheter Tip Malposition or Retraction)
- 胸管損傷 (Thoracic Duct Injury)
- 血栓 (Thromboembolism)
- 血管栓塞 (Vascular Thrombosis)
- 血管糜爛 (Vessel Erosion)

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產品型號	產品描述	Suture Plugs	Kit Type
1606052	PowerPort™ ClearVUE isp, With Attachable 6F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Microintroducer
1608052	PowerPort™ ClearVUE isp, With Attachable 8F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Microintroducer
1606062	PowerPort™ ClearVUE isp, With Attachable 6F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Intermediate
1608062	PowerPort™ ClearVUE isp, With Attachable 8F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Intermediate
1666362	PowerPort™ ClearVUE isp, With Smooth Septum and Attachable 6F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Intermediate
1668362	PowerPort™ ClearVUE isp, With Smooth Septum and Attachable 8F Polyurethane Open-Ended Single Lumen Venous Catheter	Silicone Overmold	Intermediate

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