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
Manufacturing Facility:
Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117 USA

REF **44420**

STERILE EO IUSA Rx only ② ⚠

16992-1B
22609

Manufacturer:
Medtronic, Inc.
Minneapolis, MN 55432 USA

Content components are CE approved

Strata® NSC Lumboperitoneal Shunt Kit
Strata NSC lumbo-peritoneal shuntsæt
Strata NSC lumboperitoneale shuntkit
Kit de dérivation lombopéritonéale Strata NSC
Strata NSC lumboperitoneales Shunt-Set
Κιτ οσφυοπεριτοναϊκής αναστόμωσης Strata NSC
Strata NSC lumboperitonealis sönt készlet
Kit per shunt lumboperitoneale Strata NSC
Strata NSC lumboperitonealt shuntset
Zestaw Strata NSC do wytwarzania przetoki
łędźwiowo otrzewnowej
Kit de derivação lumboperitoneal Strata NSC
Kit de derivación lumboperitoneal Strata NSC
Strata NSC lumboperitonealt shuntkit
Strata NSC Lumboperitoneal Şant Kiti

Manufacturing Facility:
Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117 USA

REF **44421**

STERILE EO IUSA Rx only ② ⚠

17137-1B
22610

Manufacturer:
Medtronic, Inc.
Minneapolis, MN 55432 USA

Content components are CE approved

REF **44430**

Strata® NSC Lumboperitoneal Shunt Kit with Integral Peritoneal Catheter
Souprava lumboperitoneálního shuntu Strata® NSC s integrovaným peritoneálním katetrem
Strata NSC lumboperitoneal shuntsæt med indbygget peritonealkateter
Strata NSC lumboperitoneale shuntkit med geintegreerde peritoneale katheter
Kit de dérivation lombopéritonéale Strata NSC avec cathéter péritonéal intégré
Strata NSC lumboperitoneales Shunt-Set mit integriertem Peritonealkatheter
Κιτ οσφυοπεριτοναϊκής αναστόμωσης Strata NSC με ολοκληρωμένο περιτοναϊκό καθετήρα
Strata NSC lumboperitoneális sönt készlet beépített peritoneális katéterrel
Kit per shunt lumboperitoneale Strata NSC con catetere peritoneale integrale
Strata NSC lumboperitonealt shuntset påmontert peritonealt kateter
Zestaw Strata NSC do wytwarzania przetoki łędźwiowo-otrzewnowej
z wbudowanym cewnikiem otrzewnowym
Kit de derivação lumboperitoneal Strata NSC com cateter peritoneal integral
Kit de derivación lumboperitoneal Strata NSC con catéter peritoneal integral
Strata NSC lumboperitonealt shuntkit med integrerad peritonealkateter
Strata NSC Lumboperitoneal Şant Kiti, Tümleşik Peritoneal Kateterli

Manufacturing Facility: Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117-5500 USA

17139-1D
24350

Manufacturer:
Medtronic, Inc., 710 Medtronic Parkway NE, Minneapolis, MN 55432-5604 USA

MR STERILE EO ② ⚠ IUSA Rx only

Content components are CE approved

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Manufacturing Facility:
 Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117 USA

REF 44465

STERILE EO !USA Rx only Ⓢ ⚠

16994-1B
 22609

Manufacturer:
 Medtronic, Inc.
 Minneapolis, MN 55432 USA

CE 0123

100. 10. 24

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105. 7. 06

Medtronic

Lumbar Catheter, Closed Tip, 84 cm
 Lumbalkateter, lukket spids, 84 cm
 Lumbale katheter, gesloten tip, 84 cm
 Cathéter lombaire, embout fermé, 84 cm
 Lumbalkatheter, geschlossene Spitze, 84 cm
 Οσφυϊκός καθετήρας, κλειστού άκρου, 84 cm
 Lumbális katéter, zárt végű, 84 cm
 Catetere lombare, a punta chiusa, da 84 cm
 Cewnik lędźwiowy z zamkniętym zakończeniem, o długości 84 cm
 Cateter lumbar, ponta fechada, 84 cm
 Catéter lumbar, punta cerrada, 84 cm
 Lumbalkateter, sluten spets, 84 cm

Manufacturing Facility:
 Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117 USA

REF 44405

STERILE EO !USA Rx only Ⓢ ⚠

16996-1B
 22609

Manufacturer:
 Medtronic, Inc.
 Minneapolis, MN 55432 USA

CE 0123

100. 10. 24

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105. 7. 06

Medtronic

Lumbar Catheter, Open Tip, 84 cm
 Lumbalkateter, åben spids, 84 cm
 Lumbale katheter, open tip, 84 cm
 Cathéter lombaire, embout ouvert, 84 cm
 Lumbalkatheter, offene Spitze, 84 cm
 Οσφυϊκός καθετήρας, 84 cm, ανοικτού άκρου
 Lumbális katéter, nyitott végű, 84 cm
 Catetere lombare, a punta aperta, da 84 cm
 Cewnik lędźwiowy z otwartą końcówką, o długości 84 cm
 Cateter lumbar, ponta aberta, 84 cm
 Catéter lumbar, punta abierta, 84 cm
 Lumbalkateter, öppen spets, 84 cm
 Lumbar Kateter, Açık Uç, 84 cm

Manufacturing Facility:
 Medtronic Neurosurgery, 125 Cremona Drive, Goleta, CA 93117 USA

REF 44406

STERILE EO !USA Rx only Ⓢ ⚠

17135-1B
 22610

Manufacturer:
 Medtronic, Inc.
 Minneapolis, MN 55432 USA

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


REF No. 44420 Strata NSC LP Shunt Kit with Closed Tip Lumbar Catheter

| |
|---|
| Included with product |
| Strata NSC Lumboperitoneal Valve |
| Lumbar Catheter, Closed Tip, Barium Impregnated |
| Peritoneal Catheter, Small Lumen, Open End, Radiopaque, with Translucent Wall, 120 cm |
| Touhy Needle, 14-Gauge, 9 cm, with Huber Tip |
| Fixation Tab, Lumbar (2) |
| Fixation Tab, Large (1) |
| Blunt Needle, 20-Gauge |
| Guidewire with Adjustable Stop |
| Strain Relief |

REF No. 44430 Strata NSC LP Shunt Kit with Closed Tip Lumbar Catheter and Integral Peritoneal Catheter, Small Lumen

| |
|--|
| Included with product |
| Strata NSC Lumboperitoneal Valve with Integral Lumbar Catheter, Closed Tip, Barium Impregnated |
| Peritoneal Catheter, Small Lumen, Open End, Radiopaque, with Translucent Wall, 120 cm |
| Touhy Needle, 14-Gauge, 9 cm, with Huber Tip |
| Fixation Tab, Lumbar (2) |
| Fixation Tab, Large (1) |
| Blunt Needle, 20-Gauge |
| Guidewire with Adjustable Stop |
| Strain Relief |

REF No. 44421 Strata NSC LP Shunt Kit with Open Tip Lumbar Catheter

| |
|---|
| Included with product |
| Strata NSC Lumboperitoneal Valve |
| Lumbar Catheter, Open Tip, Barium Impregnated |
| Peritoneal Catheter, Small Lumen, Open End, Radiopaque, with Translucent Wall, 120 cm |
| Touhy Needle, 14-Gauge, 9 cm, with Huber Tip |
| Fixation Tab, Lumbar (2) |
| Fixation Tab, Large (1) |
| Strain Relief |

REF No. 44465 Strata NSC LP Valve

| |
|-----------------------|
| Included with product |
| Strain Relief |

REF No. 44405 Lumbar Catheter Kit, Closed Tip, 80 cm

| |
|---|
| Included with product |
| Lumbar Catheter, Closed Tip, Barium Impregnated |
| Touhy Needle, 14-Gauge, 9 cm, with Huber Tip |
| Fixation Tab, Lumbar (2) |
| Blunt Needle, 20-Gauge |
| Guidewire with Adjustable Stop |
| Strain Relief |

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REF No. 44406 Lumbar Catheter Kit, Open Tip, 80 cm

| |
|---|
| Included with product |
| Lumbar Catheter, Open Tip, Barium Impregnated |
| Touhy Needle, 14-Gauge, 9 cm, with Huber Tip |
| Fixation Tab, Lumbar (2) |
| Blunt Needle, 20-Gauge |
| Strain Relief |

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REF No. 22063 Catheter Connector, Stepdown

| |
|--|
| Included with product |
| Silicone Tubing, Standard, 5cm, Barium Impregnated |
| 此導管連接管適用於連接 PS Medical 內徑 0.7mm 及 1.3mm 的導管。 |

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