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“AMO” TECNIS Symphony Extended Range of Vision 1-Piece IOL

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
Johnson & Johnson VISION

“TECNIS” prints 4-color
Process of Symfony Gold

TECNIS
Symfony®
Extended Range
of Vision IOL

TECNIS and TECNIS Symfony are trademarks of
Johnson & Johnson Surgical Vision, Inc.
Patents: www.jiv-patents.com
© Johnson & Johnson Surgical Vision, Inc. 2018

TECNIS
Symfony®
Extended Range
of Vision IOL

TECNIS Symfony® Extended Range of Vision IOL

Johnson & Johnson VISION

TECNIS
Symfony®
Extended Range
of Vision IOL

STERILE EO

Etienne Oxide

Do Not Reuse

Do Not Reseal

Do Not Sterilize

Do Not Heat

Do Not Expose to Moisture

Do Not Use if Damaged

Keep Away From Sunlight

Keep Ambient Temperature

Upper Limit 45°C (113°F)

Lower Limit 5°C (41°F)

Consist Instructions For Use

AMO Ireland
Block B
Lilly Valley Office Campus
Quarway, Co. Dublin, Ireland

TECNIS
Symfony®
Extended Range
of Vision IOL

EN The package contains one (1) Dimpled Aspheric Extended Range of Vision IOL. Data regarding refractive outcomes for this IOL are available at www.jiv.com.

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製造廠名: (O)Johnson & Johnson Surgical Vision, Inc. (P)AMO Puerto Rico Manufacturing, Inc.


製造廠址: (O)1700 East St. Andrew Place Santa Ana, CA USA 92705
(P)Road 402 North, Anasco, Puerto Rico 00610, U.S.A.

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Figure 1
Light Transmittance

LEGEND:

Spectral transmittance curve of a typical 5-diopter IOL (thinnest), UV cut-off at 10%T is 374 nm
Spectral transmittance curve of a typical 20-diopter IOL, UV cut-off at 10%T is 376 nm
Spectral transmittance curve of a typical 34-diopter IOL (thickest), UV cut-off at 10%T is 375 nm
Spectral transmittance curve of a 53-year-old phakic eye, from Boelker, E.A., and Waller J.R. *Transmission of the Ocular Media. Investigative Ophthalm.* 1962;1:776-783.

Note: The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance of the IOLs (5-34D) made with material. Spectral transmission measurements were taken in water at room temperature (ref: TR7476).

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